

The Ultimate in Performance Veterinary Power Tools





VDRIVE™ V-MBU SMALL BONE POWERED INSTRUMENT RANGE



V-MBU-470

Fitted with sagittal saw attachment and A.O. drill.

Advanced function controls

- High speed and oscillating drilling.
- Low speed drilling, screw driving and tapping.
- Torque limiting for screw driving.
- Forward and reverse triggers.
- Mode selector button with mechanical safety lock.

Power options

- 14.6V Li-lon aseptic batteries.
- 13.2V Li-lon sterile batteries.
- Corded direct power.

V-MBU SYSTEM

HANDPIECE




| ORDER CODE | MODEL | DESCRIPTION | SPEED (rpm) |
|------------|-----------|--------------|----------------|
| 1292444 | V-MBU-470 | Twin Trigger | 0 - 1350 |

ATTACHMENTS

| ORDER CODE | MODEL | DESCRIPTION | SPEED (rpm) | CANNULATION Ø mm |
|--------------------------|----------|--|--|---------------------|
| 18780 | V-DU-470 | 0.5 - 4.0 mm (5/32") Keyed | 0 - 1350 | 4.1 |
| 18770 | V-DU-470 | 0.5 - 7.4 mm (1/4") Keyed | 0 - 1350 | 4.2 |
| 18800 | V-DU-470 | A.O. Synthes/ASIF | 0 - 1350 | 2.55 |
| 18790 | V-DU-470 | Hudson/Zimmer | 0 - 1350 | 4.2 |
| 19480 | V-BU-472 | High speed bur Uses R Series burs | 0 - 27000 | - |
| Standard b Order code | | Long bur guard Order code 18730 | Extra long bur guard Order code 18740 | |
| 18840 | V-RU-470 | Hudson/Zimmer | 0 - 350 | 4.2 |
| 18830 | V-RU-470 | A.O. Synthes | 0 - 350 | 4.2 |
| 18820 | V-RU-470 | 0.5 - 7.4 mm (1/4") | 0 - 350 | 4.2 |
| 18870 | V-SU-470 | Mini Sagittal saw Uses S85 type blades | 0 - 22000 cpm | - |
| 18880 | V-NU-470 | Quick release sagittal saw Uses S88 type blades | 0 - 20500 cpm | - |
| 18890 | V-CU-470 | Small reciprocating saw Uses S92 type blades | 0 - 8000 cpm | - |







VDRIVE™ V-MBQ LARGE BONE POWERED INSTRUMENT RANGE



HANDDIECE SPECIFICATION

| HANDPIECE SPECIFI | CATION | | | |
|-------------------|----------------|---------------------|---------------------------|----------------|
| | ORDER CODE | MODEL | DESCRIPTION | SPEED (rpm) |
| Valrice | 1291774 | V-MBQ-707 | Single Trigger | 0 - 1350 |
| Valence | 1291884 | V-MBQ-708 | Twin Trigger | 0 - 1350 |
| | * only to be u | sed with De Soutter | Medical V-MBQ attachment: | s |



| | ORDER CODE | MODEL | DESCRIPTION | SPEED (cpm) |
|-------|------------|-----------|--|----------------|
| lifes | 1293884 | V-DBO-800 | Dedicated straight case TPLO saw | 0 - 13500 |
| | | | Uses De Soutter Medical R series b * Also accepts Synthes® blades | lade |

MODULAR HANDPIECE SPECIFICATION

| | ORDER CODE | MODEL | DESCRIPTION | SPEED (rpm) |
|---|------------|------------|--|----------------------------------|
| 0 | 1293774 | V-MBQ-807* | Modular straight case handpiece * To be used with De Soutter Medic | 0 - 1350 al V-MBQ attachments |

V-MBQ ATTACHMENTS

| | ORDER CODE | MODEL | DESCRIPTION | SPEED (rpm) | CANNULATION Ø MM |
|-----|------------|-----------|--|------------------------|---------------------|
| | 18320 | V-DQ-708 | 0.5 - 6.4 mm (1/4") keyed | 0 - 1350 | 4.4 |
| | 18330 | V-DQ-708 | Small A.O Synthes/ASIF | 0 - 1350 | 2.55 |
| | 18440 | V-RQ-708 | 0.5 - 6.4 mm (1/4") keyed | 0 - 330 | 4.4 |
| | 18290 | V-KQ-707 | Sagittal saw Uses S89 type blades | 0 - 13500 cpm | ٦ |
| | 18300 | V-NQ-707 | Small sagittal saw Uses S88 type blades | 0 - 20000 cpr | m |
| 120 | 18690 | V-0Q-708* | TPLO oscillating saw Uses De Soutter Medical R series k * Also accepts Synthes® blades | 0 - 13500 cpm blade | ١ |
| | 18270 | V-0Q-708 | TPLO oscillating saw Uses Slocum blades | 0 - 13500 cpm | ٦ |
| | 18310 | V-WQ-707 | 0.7 - 4.0 mm wire capacity | 0 - 1350 | 4.0 |

BATTERY AND ELECTRICAL COMPONENTS

Why choose Lithium Ion battery technology?

De Soutter Medical Lithium Ion batteries offer the surgical team and sterile services a number of tangible benefits over conventional battery technology, including:

- Li-lon technology improves battery 'capacity to weight' ratio by around 100% when compared to NiMH or NiCad chemistries.
- Significantly extends the usable run time of the power tool.
- Eliminates 'memory effect' normally associated with other battery types.
- Unique **fuel gauge** incorporated within De Soutter Medical Li-Ion batteries allows the surgical team to instantly assess the battery capacity
 - both before and during use.



STERILE BATTERIES - LI-ION

| ORDER CODE | MODEL | DESCRIPTION | VOLTAGE | CAPACITY (mAh) |
|------------|--------|---------------|---------|----------------|
| 17210 | SB-703 | Small battery | 13.2V | 1100 |
| 17220 | SB-704 | Large battery | 13.2V | 2500 |

ASEPTIC BATTERIES



SMALL BATTERY

Battery pack



Battery housing



Battery shield

| ORDER CODE | MODEL | DESCRIPTION | VOLTAGE | CAPACITY (mAh) |
|------------|----------|-----------------|---------|----------------|
| 17910 | AB-703 | Battery pack | 10.95V | 2000 |
| 18850 | V-AH-703 | Battery housing | - | - |
| 18100 | AS-703 | Battery shield | - | - |

LARGE BATTERY



Battery pack



Battery housing



Battery shield

| ORDER CODE | MODEL | DESCRIPTION | VOLTAGE | CAPACITY (mAh) |
|------------|----------|-----------------|---------|----------------|
| 17920 | AB-704 | Battery pack | 14.6V | 2000 |
| 18860 | V-AH-704 | Battery housing | - | - |
| 18110 | AS-704 | Battery shield | - | - |

BATTERY AND ELECTRICAL COMPONENTS

BC-700 MULTI BAY CHARGER AND MODULE OPTIONS

| CHARGER ORDER CODE | CHARGER MODULE DESCRIPTION | CHARGER MODULE ORDER CODE | FOR BATTERY TYPE |
|-----------------------|-----------------------------------|------------------------------|------------------|
| 15700 | Without charger module fitted | | |
| 16230 | With 4-off CA-703 charger modules | 18120 | AB-703 / AB-704 |
| 16220 | With 4-off CA-700 charger modules | 15750 | SB-703 / SB-704 |

SINGLE BAY CHARGERS

| | ORDER CODE | MODEL | FOR BATTERY TYPE |
|------|------------|--------|------------------|
| -(> | 17990 | BC-706 | SB-703/SB-704 |
| | 18000 | BC-707 | AB-703 |
| | 18010 | BC-708 | AB-704 |

MAINS CABLES

Mains power, 3 metres long

| ORDER CODE | COUNTRY |
|------------|----------------------------|
| 8290 | UK |
| 8300 | Euro |
| 8310 | US / Canada |
| 8320 | Australia / New Zealand |
| 14700 | Swiss |

BC-700 BATTERY CHARGER

The BC-700 Battery Charger requires minimal user input, simply insert battery and the BC-700 will charge, recover and recondition it to its optimum state.



TPLO SAW BLADES

De Soutter Medical manufactured TPLO saw blades offer a range of advantages;

- Thinner blades reduces weight, vibration and requires less cutting force
- Reduced cut thickness minimises heat gain
- A sharper blade lasts longer
- Tooth offset improves cutting performance
- Tapered hub design reduces vibration and eliminates risk of blade coming undone
- Simple single captive screw makes fitting easier



Synthes® compatible

| | в |
|--------------|----|
| | ν. |
| | |
| | |
| | |
| | |
| The state of | |
| 1000 | |
| | |
| | |
| | |
| 7007 | |
| 100 | |
| | |
| | |

| | Length 45mm | | | | Length | 50mm | |
|--------------|-------------|-------|-------|-------|--------|-------|-------|
| Radius | 12 | 15 | 18 | 21 | 24 | 27 | 30 |
| DSM/Synthes® | 20160 | 20170 | 20180 | 20190 | 20200 | 20210 | 20220 |
| Slocum | 20070 | 20080 | 20090 | 20100 | 20110 | 20120 | 20130 |

Synthes® is a registered trademark of Synthes-Stratec. De Soutter Medical are not associated with the Synthes® trademark.

STERILISATION

| | SYSTEM | Half size Container & insert tray Order code | Half size Insert tray Order code |
|---|--|---|---|
| | V-MBU ASEPTIC | 19210 | 19190 |
| | V-MBQ-70X ASEPTIC | 19230 | 19220 |
| | V-MBQ-807 ASEPTIC | 20430 | 20380 |
| | UNIVERSAL | 19000 | 667923 |
| Sterile filtered container Accepts half size wire insert tray Order code: 19370 | Filter paper - single use pack of 100 Order code: 657243 | Security tag, pack of 100 Order code: 657253 | Sterile indicator label, pack of 100 Order code: 657263 Lid gasket Order code: 669903 |

CORDLESS BATTERY POWERED SAW

SAW SPECIFICATION



| ORDER CODE | MODEL | VOLTAGE | OSCILLATING SPEED (CPM) | WEIGHT (GRAMS) |
|------------|---------|----------|----------------------------|-------------------|
| 1286224 | CSB-100 | 10.8 Vdc | 16000 | 1050 |

Supplied with saw blade part no. 6250.

BATTERY AND BATTERY CHARGER

| ORDER CODE | MODEL | VOLTAGE | SUPPLY FREQUENCY |
|------------|--------|---------------|------------------|
| 15020 | BC-100 | 100 - 240 Vac | 50 - 60 Hz |
| 15010 | | 10.8 Vdc NiMH | |

CSB-100-A CAST SAW SYSTEM



| DESCRIPTION | ORDER CODE |
|--------------------|------------|
| CSB-100 Saw | 1286224 |
| CB-100 Battery x 2 | 15010 |
| BC-100 Charger | 15020 |
| Saw blade | 6250 |
| Hexagon wrench | 204903 |
| Carry case | 15030 |

CSB-100

The De Soutter Medical CSB-100 saw allows the user total freedom of movement, which is particularly useful for the veterinary profession where the cordless concept is the only practical and affordable solution.





United Kingdom

(International Head Office) De Soutter Medical Limited Halton Brook Business Park, Weston Road, Aston Clinton, Aylesbury, Bucks, HP22 5WF Tel: +44 (0) 1296 634 000 Fax: +44 (0) 1296 634 033

Email: info@de-soutter.com

Belgium

De Soutter Medical Belgium Bessemerstraat 14, 3620 Lanaken Tel: +32 (0) 89/47 15 37 Fax: +32 (0) 89/70 12 19 Email: belgium@de-soutter.com

Italy

De Soutter Medical Italia Località Forna ce SNC 27022 Casorate Primo - PV Tel: +39 (0) 2 9009 4098 Fax: +39 (0) 2 9009 2673 Email: italy@de-soutter.com

Australia

De Soutter Medical Australia Pty Ltd 2/12-14 Apollo Drive, Hallam, Victoria 3803 Tel: +61 (0) 3 9702 4441

Fax: +61 (0) 3 9702 4484 Email: australia@de-soutter.com

France

De Soutter Medical France 1 252 Avenue Parc des Expositions 33260 La Teste de Buch Tel: +33 (0) 5 56 54 89 36 Fax: +33 (0) 9 70 61 37 60 Email: france@de-soutter.com

Netherlands

De Soutter Medical Netherlands Gelderlandhaven 2X 3433 PG Nieuwegein Tel: +31 (0) 85 0491480 Fax: +31 (0) 85 0491489 Email: nederland@de-soutter.com

Austria

De Soutter Medical Austria Zweigniederlassung Österreich Dietrichsteingasse 10, A-3400 Klosterneuburg Tel: +43 (0) 676 96 71 770 Fax: +43 (0) 2243 21 656 Email: austria@de-soutter.com

211-EN-11

Germany

De Soutter Medical Germany Bahnhofstraße 4, Kladenfloss D-66625 Nohfelden Tel: +49 (0) 68 52-99 12 46 Fax: +49 (0) 68 52-99 12 47 Email: deutschland@de-soutter.com

United States of America

De Soutter Medical USA, Inc. 224 Rolling Hill Road, Suite 12A Mooresville, NC 28117 Tel: +1 (704) 655 9040 Fax: +1 (704) 987 2035 Email: usa@de-soutter.com

